

KIAA0391 Polyclonal Antibody

Catalog Number:E-AB-62319



Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Rat
Immunogen	Recombinant fusion protein of human KIAA0391 (NP_055487.2).
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Formulation	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Applications Recommended Dilution

WB	1:500-1:2000
-----------	--------------

Data



Western blot analysis of extracts of various cell lines using KIAA0391 Polyclonal Antibody at dilution of 1:1000.

Observed Mw:67kDa

Calculated Mw:24kDa/56kDa/65kDa/67kDa

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

PRORP (Protein Only RNase P Catalytic Subunit) is a Protein Coding gene. Diseases associated with PRORP include Perrault Syndrome 1 and Syndromic X-Linked Intellectual Disability Type 10. Among its related pathways are tRNA processing and Gene Expression. An important paralog of this gene is ENSG00000258790.

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017