

THAP2 Polyclonal Antibody

Catalog No. E-AB-65443

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Mouse,Rat
Immunogen	Recombinant fusion protein of human THAP2 (NP_113623.1).
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Applications Recommended Dilution

WB 1:500-1:2000

Data



Western blot analysis of extracts of Mouse small intestine using THAP2 Polyclonal Antibody at dilution of 1:3000.

Observed Mw:30kDa
Calculated Mw:26kDa

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

THAP2 (THAP Domain Containing 2) is a Protein Coding gene. Gene Ontology (GO) annotations related to this gene include nucleic acid binding. An important paralog of this gene is ENSG00000258064.

For Research Use Only