

CAPZA1 Polyclonal Antibody

Catalog No. E-AB-63755

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human,Mouse,Rat
Immunogen	Recombinant protein of human CAPZA1
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.

Applications Recommended Dilution

WB	1:200-1:1000
-----------	--------------

Data



Western blot analysis of extracts of various cell lines
with CAPZA1 Polyclonal Antibody

Observed Mw:38kDa
Calculated Mw:32kDa

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

CAPZA1 is a member of the F-actin capping protein alpha subunit family. This gene encodes the alpha subunit of the barbed-end actin binding protein. The protein regulates growth of the actin filament by capping the barbed end of growing actin filaments.

For Research Use Only