

OSBP Polyclonal Antibody

Catalog No. E-AB-63770

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Description


Reactivity	Human, Mouse, Rat
Immunogen	Recombinant protein of human OSBP
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.

Applications

Recommended Dilution

WB 1:500-1:2000

Data



Western blot analysis of extracts of various cell lines
with OSBP Polyclonal Antibody

Observed Mw:110kDa

Calculated Mw:89kDa

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

Oxysterol binding protein is an intracellular protein that is believed to transport sterols from lysosomes to the nucleus where the sterol down-regulates the genes for the LDL receptor, HMG-CoA reductase, and HMG synthetase

For Research Use Only