

# LARP1 Polyclonal Antibody

Catalog Number:E-AB-65483



**Note:** Centrifuge before opening to ensure complete recovery of vial contents.


## Description

<b>Reactivity</b>	Human
<b>Immunogen</b>	Recombinant Protein of human LARP1
<b>Host</b>	Rabbit
<b>Isotype</b>	IgG
<b>Purification</b>	Affinity purification
<b>Conjugation</b>	Unconjugated
<b>Formulation</b>	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Applications Recommended Dilution

<b>WB</b>	1:500-1:2000
-----------	--------------

## Data



Western blot analysis of extracts of various cell lines using LARP1 Polyclonal Antibody at dilution of 1:1000.

**Observed Mw:163kDa**  
**Calculated Mw:116kDa/123kDa**

## Preparation & Storage

**Storage** Store at -20°C. Avoid freeze / thaw cycles.

## Background

LARP1 (La Ribonucleoprotein 1, Translational Regulator) is a Protein Coding gene. Diseases associated with LARP1 include Hepatocellular Carcinoma. Among its related pathways are Translational Control. Gene Ontology (GO) annotations related to this gene include translation initiation factor binding. An important paralog of this gene is LARP1B.

## For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: [www.elabscience.com](http://www.elabscience.com)

Tel: 1-832-243-6086

Email: [techsupport@elabscience.com](mailto:techsupport@elabscience.com)

Fax: 1-832-243-6017