

LARP1 Polyclonal Antibody

Catalog Number:E-AB-65483



Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human
Immunogen	Recombinant Protein of human LARP1
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Formulation	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Applications Recommended Dilution

WB	1:500-1:2000
-----------	--------------

Data



Western blot analysis of extracts of various cell lines using LARP1 Polyclonal Antibody at dilution of 1:1000.

Observed Mw:163kDa

Calculated Mw:116kDa/123kDa

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

LARP1 (La Ribonucleoprotein 1, Translational Regulator) is a Protein Coding gene. Diseases associated with LARP1 include Hepatocellular Carcinoma. Among its related pathways are Translational Control. Gene Ontology (GO) annotations related to this gene include translation initiation factor binding. An important paralog of this gene is LARP1B.

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017