

UQCR10 Polyclonal Antibody

Catalog No. E-AB-62432

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human,Mouse,Rat
Immunogen	Recombinant fusion protein of human UQCR10
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide,50% glycerol,pH7.3.

Applications Recommended Dilution

WB 1:500-1:2000

Data



Western blot analysis of extracts of various cell lines using UQCR10 Polyclonal Antibody at 1:1000 dilution.

Observed Mw:7kDa

Calculated Mw:6kDa/7kDa

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

UCRC is a subunit of mitochondrial complex III (ubiquinol-cytochrome c reductase; EC 1.10.2.2), which forms the middle segment of the respiratory chain of the inner mitochondrial membrane (Schagger et al., 1995 [PubMed 8592474]).[supplied by OMIM, Mar 2008]

For Research Use Only