

TSR1 Polyclonal Antibody

Catalog No. E-AB-65500

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human, Mouse
Immunogen	Recombinant fusion protein of human TSR1 (NP_060598.3).
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Applications Recommended Dilution

WB 1:500-1:2000

Data



Western blot analysis of extracts of various cell lines using TSR1 Polyclonal Antibody at dilution of 1:3000.

Observed Mw:110kDa

Calculated Mw:91kDa

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

TSR1 (TSR1 Ribosome Maturation Factor) is a Protein Coding gene. Diseases associated with TSR1 include Epithelial Malignant Thymoma and Dendritic Cell Thymoma. Among its related pathways are rRNA processing in the nucleus and cytosol and Gene Expression. An important paralog of this gene is BMS1.

For Research Use Only