

MT-ND1 Polyclonal Antibody

Catalog No. E-AB-62128

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human, Mouse, Rat
Immunogen	A synthetic peptide of mouse MT-ND1
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.

Applications Recommended Dilution


WB 1:500 - 1:1000

Data



Western blot analysis of extracts of various cell lines with MT-ND1 Polyclonal Antibody

Observed Mw: 26-38kDa
Calculated Mw: 36kDa



Western blot analysis of extracts of various cell lines with MT-ND1 Polyclonal Antibody

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

For Research Use Only