# MPG Polyclonal Antibody

Catalog Number: E-AB-60952



*Note:* Centrifuge before opening to ensure complete recovery of vial contents.

Description	
Reactivity	Human,Mouse
Immunogen	Recombinant fusion protein of human MPG (NP_001015052.1).
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Formulation	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Applications	Recommended Dilution
WB	1:500-1:2000
IF	1:50-1:100
Data	



Western blot analysis of extracts of various cell lines using MPG Polyclonal Antibody at dilution of 1:1000.

Observed Mw:37kDa Calculated Mw:30kDa/32kDa



Immunofluorescence analysis of GFP-RNF168 transgenic U2OS cells using MPG Polyclonal Antibody

Immunofluorescence analysis of U2OS cells using MPG Polyclonal Antibody



Immunofluorescence analysis of A549 cells using MPG Polyclonal Antibody

## **Preparation & Storage**

#### Storage

Store at -20°C. Avoid freeze / thaw cycles.

#### **Background**

Maintenance of DNA sequences is necessary for vertebrates and other life. DNA is under constant stress by a plethora of DNA-damaging agents present in both the environment and within cells. The potentially deleterious effects of DNA

### For Research Use Only

A Reliable Research Partner in Life Science and Medicine

## MPG Polyclonal Antibody

Catalog Number: E-AB-60952



lesions in cells are elegantly resolved by sophisticated DNA repair systems, including base excision repair (BER), nucleotide excision repair (NER) and DNA repair methyltransferase (MTase). Methylated bases, such as 3-methyladenine (3MeA) and 7-methylguanine (7MeG) can be formed by agents in the environment and by endogenous cellular processes. Consequently, in the absence of exposure to environmental agents, DNA methylation damage can be incurred on the genomic DNA of normal mammalian cells. DNA N-glycosylases are base excision-repair proteins that locate and cleave damaged bases from DNA as the first step in restoring the sequence.

### For Research Use Only

A Reliable Research Partner in Life Science and Medicine
Toll-free: 1-888-852-8623 Tel: 1-832-243-6086
Web: www.elabscience.com Email: techsupport@elabscience.com