

GGH Polyclonal Antibody

Catalog No. E-AB-62983

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Description


Reactivity	Human,Mouse
Immunogen	Recombinant fusion protein of human GGH
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide,50% glycerol,pH7.3.

Applications

Recommended Dilution

WB 1:500-1:2000

Data



Western blot analysis of extracts of various cell lines using GGH Polyclonal Antibody at 1:1000 dilution.

Observed Mw:32kDa
Calculated Mw:35kDa

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

This gene catalyzes the hydrolysis of foylpoly-gamma-glutamates and antifoylpoly-gamma-glutamates by the removal of gamma-linked polyglutamates and glutamate.

For Research Use Only