

## KLK7 Polyclonal Antibody

**Catalog No.** E-AB-61000

**Note:** Centrifuge before opening to ensure complete recovery of vial contents.


### Description

|                     |                                                      |
|---------------------|------------------------------------------------------|
| <b>Reactivity</b>   | Human, Mouse                                         |
| <b>Immunogen</b>    | Recombinant protein of human KLK7                    |
| <b>Host</b>         | Rabbit                                               |
| <b>Isotype</b>      | IgG                                                  |
| <b>Purification</b> | Affinity purification                                |
| <b>Conjugation</b>  | Unconjugated                                         |
| <b>Buffer</b>       | PBS with 0.02% sodium azide and 50% glycerol pH 7.4. |

### Applications Recommended Dilution

**WB** 1:500-1:2000

### Data



Western blot analysis of extracts of various cell lines  
with KLK7 Polyclonal Antibody

**Observed Mw:30kDa**

**Calculated Mw:19kDa/27kDa**

### Preparation & Storage

**Storage** Store at -20°C. Avoid freeze / thaw cycles.

### Background

This gene encodes a member of the kallikrein subfamily of serine proteases. These enzymes have diverse physiological functions and many kallikrein genes are biomarkers for cancer. The encoded protein has chymotrypsin-like activity and plays a role in the proteolysis of intercellular cohesive structures that precedes desquamation, the shedding of the outermost layer of the epidermis. The encoded protein may play a role in cancer invasion and metastasis, and increased expression of this gene is associated with unfavorable prognosis and progression of several types of cancer. Polymorphisms in this gene may play a role in the development of atopic dermatitis. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene, which is one of fifteen kallikrein subfamily members located in a gene cluster on chromosome 19.

### For Research Use Only