

SMYD2 Polyclonal Antibody

Catalog No. E-AB-61990

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human,Mouse
Immunogen	Recombinant fusion protein of human SMYD2 (NP_064582.2).
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Applications Recommended Dilution

IF 1:10-1:100

Data



Immunofluorescence analysis of U2OS cells using
SMYD2 Polyclonal Antibody

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

SMYD2 (SET And MYND Domain Containing 2) is a Protein Coding gene. Among its related pathways are Lysine degradation and PKMTs methylate histone lysines. Gene Ontology (GO) annotations related to this gene include p53 binding and protein-lysine N-methyltransferase activity. An important paralog of this gene is SMYD1.

For Research Use Only