

BCL2L14 Polyclonal Antibody

Catalog No. E-AB-61992

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human, Mouse, Rat
Immunogen	Recombinant fusion protein of human BCL2L14 (NP_620049.1).
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Applications Recommended Dilution


WB 1:500-1:2000 IF
1:10-1:100

Data



Western blot analysis of extracts of various cell lines using BCL2L14 Polyclonal Antibody at dilution of 1:1000.

Observed Mw:39kDa
Calculated Mw:28kDa/30kDa/36kDa



Immunofluorescence analysis of HeLa cells using BCL2L14 Polyclonal Antibody

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

The protein encoded by this gene belongs to the BCL2 protein family. BCL2 family members form hetero- or homodimers and act as anti- or pro-apoptotic regulators that are involved in a wide variety of cellular activities. Overexpression of this gene has been shown to induce apoptosis in cells. Three alternatively spliced transcript variants encoding two distinct isoforms have been reported for this gene.

For Research Use Only