

HAGH Polyclonal Antibody

Catalog No. E-AB-61250

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human,Mouse,Rat
Immunogen	Recombinant fusion protein of human HAGH
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide,50% glycerol,pH7.3.

Applications Recommended Dilution


WB	1:500-1:2000
IF	1:10-1:100

Data



Western blot analysis of extracts of various cell lines using HAGH Polyclonal Antibody at 1:1000 dilution.

Observed Mw:34kDa
Calculated Mw:25kDa/28kDa/33kDa



Immunofluorescence analysis of A549 cells using HAGH Polyclonal Antibody Blue: DAPI for nuclear staining.

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

The enzyme encoded by this gene is classified as a thiolesterase and is responsible for the hydrolysis of S-lactoyl-glutathione to reduced glutathione and D-lactate. Three transcript variants encoding different isoforms have been found for this gene.

For Research Use Only