

KHDRBS3 Polyclonal Antibody

Catalog No. E-AB-63045

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human, Mouse
Immunogen	Recombinant fusion protein of human KHDRBS3 (NP_006549.1).
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Applications Recommended Dilution

IF 1:50-1:100

Data



Immunofluorescence analysis of MCF7 cells using KHDRBS3 Polyclonal Antibody

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

KHDRBS3 (KH RNA Binding Domain Containing, Signal Transduction Associated 3) is a Protein Coding gene. Diseases associated with KHDRBS3 include Extratemporal Epilepsy and Subacute Leukemia. Among its related pathways are Signaling by GPCR and Signaling by PTK6. Gene Ontology (GO) annotations related to this gene include nucleic acid binding and RNA binding. An important paralog of this gene is KHDRBS2.

For Research Use Only