

TPPP3 Polyclonal Antibody

Catalog No. E-AB-61356

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human, Mouse, Rat
Immunogen	Recombinant fusion protein of human TPPP3 (NP_057224.2).
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Applications Recommended Dilution


WB 1:500-1:2000 IHC
1:50-1:200

Data




Western blot analysis of extracts of various cell lines using TPPP3 Polyclonal Antibody at dilution of 1:1000.

Observed Mw:19kDa
Calculated Mw:18kDa



Immunohistochemistry of paraffin-embedded Human gastric cancer using TPPP3 Polyclonal Antibody at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded Human colon carcinoma using TPPP3 Polyclonal Antibody at dilution of 1:100 (40x lens).

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

For Research Use Only

Background

TPPP3 (Tubulin Polymerization Promoting Protein Family Member 3) is a Protein Coding gene. Diseases associated with TPPP3 include Cerebrospinal Fluid Leak and Creutzfeldt-Jakob Disease. Gene Ontology (GO) annotations related to this gene include tubulin binding. An important paralog of this gene is TPPP.

For Research Use Only