

POLR2H Polyclonal Antibody

Catalog No. E-AB-61423

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human, Mouse, Rat
Immunogen	Recombinant fusion protein of human POLR2H (NP_006223.2).
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Applications Recommended Dilution

WB 1:500-1:2000

Data



Western blot analysis of extracts of various cell lines using POLR2H Polyclonal Antibody at dilution of 1:1000.

Observed Mw:17kDa

Calculated Mw:9kDa/12kDa/13kDa/17kDa/19kDa

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

The three eukaryotic RNA polymerases are complex multisubunit enzymes that play a central role in the transcription of nuclear genes. This gene encodes an essential and highly conserved subunit of RNA polymerase II that is shared by the other two eukaryotic DNA-directed RNA polymerases, I and III. Alternative splicing results in multiple transcript variants of this gene.

For Research Use Only