

C16orf80 Polyclonal Antibody

Catalog Number:E-AB-61522

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human,Mouse,Rat
Immunogen	Recombinant fusion protein of human C16orf80 (NP_037374.1).
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Formulation	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Applications Recommended Dilution


WB	1:500-1:2000
IF	1:50-1:100

Data



Western blot analysis of extracts of various cell lines using C16orf80 Polyclonal Antibody at dilution of 1:1000.

Observed Mw:23kDa
Calculated Mw:22kDa



Immunofluorescence analysis of A549 cells using C16orf80 Polyclonal Antibody

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

CFAP20 (Cilia And Flagella Associated Protein 20) is a Protein Coding gene. Diseases associated with CFAP20 include Familial Atrial Fibrillation. An important paralog of this gene is CFAP20DC.

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017