

TXNL4B Polyclonal Antibody

Catalog No. E-AB-61535

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human, Mouse
Immunogen	Recombinant fusion protein of human TXNL4B (NP_060323.1).
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Applications Recommended Dilution


WB 1:500-1:2000 IF
1:50-1:100

Data



Western blot analysis of extracts of Mouse liver using TXNL4B Polyclonal Antibody at dilution of 1:1000.

Observed Mw:17kDa
Calculated Mw:17kDa



Immunofluorescence analysis of A-549 cells using TXNL4B Polyclonal Antibody

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

TXNL4B (Thioredoxin Like 4B) is a Protein Coding gene. Diseases associated with TXNL4B include Anhaptoglobinemia and Epithelial-Stromal Tgfbi Dystrophy. An important paralog of this gene is TXNL4A.

For Research Use Only