

TRMT1 Polyclonal Antibody

Catalog No. E-AB-62482

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human, Mouse
Immunogen	Recombinant fusion protein of human TRMT1
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Applications Recommended Dilution

WB 1:500-1:2000

Data



Western blot analysis of extracts of various cell lines using TRMT1 Polyclonal Antibody at 1:1000 dilution.

Observed Mw:90kDa

Calculated Mw:69kDa/72kDa

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

Dimethylates a single guanine residue at position 26 of most tRNAs using S-adenosyl-L-methionine as donor of the methyl groups.

For Research Use Only