

KLF9 Polyclonal Antibody

Catalog No. E-AB-61578

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human,Mouse,Rat
Immunogen	Recombinant fusion protein of human KLF9
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.01% thiomersal,50% glycerol,pH7.3.

Applications Recommended Dilution


WB 1:500-1:2000

Data



Western blot analysis of extracts of various cell lines using KLF9 Polyclonal Antibody at 1:1000 dilution.

Observed Mw:30KDa
Calculated Mw:27kDa



Western blot analysis of extracts of various cell lines using KLF9 Polyclonal Antibody at 1:1000 dilution.

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

The protein encoded by this gene is a transcription factor that binds to GC box elements located in the promoter. Binding of the encoded protein to a single GC box inhibits mRNA expression while binding to tandemly repeated GC box elements activates transcription.

For Research Use Only