

NSUN6 Polyclonal Antibody

Catalog No. E-AB-61583

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Description


Reactivity	Human, Mouse
Immunogen	Recombinant protein of human NSUN6
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.

Applications Recommended Dilution

WB 1:500 - 1:2000 IP


1:50 - 1:200

Data



Western blot analysis of extracts of various cell lines with NSUN6 Polyclonal Antibody

Observed Mw: 52kDa
Calculated Mw: 51kDa



Immunoprecipitation analysis of 200ug extracts of HeLa cells with NSUN6 Polyclonal Antibody

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

For Research Use Only