

PHYHD1 Polyclonal Antibody

Catalog No. E-AB-61585

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human,Mouse,Rat
Immunogen	Recombinant fusion protein of human PHYHD1 (NP_001094346.1).
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Applications Recommended Dilution

IF	1:50-1:100
-----------	------------

Data



Immunofluorescence analysis of A-549 cells using
PHYHD1 Polyclonal Antibody

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

PHYHD1 (Phytanoyl-CoA Dioxygenase Domain Containing 1) is a Protein Coding gene. Gene Ontology (GO) annotations related to this gene include dioxygenase activity.

For Research Use Only