

NGFRAP1 Polyclonal Antibody

Catalog No. E-AB-61628

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human, Mouse
Immunogen	Recombinant fusion protein of human NGFRAP1 (NP_996798.1).
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Applications Recommended Dilution

IF 1:50-1:100

Data



Immunofluorescence analysis of MCF-7 cells using
NGFRAP1 Polyclonal Antibody

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

Protein BEX3 is a protein that in humans is encoded by the NGFRAP1 gene. May be a signaling adapter molecule involved in p75NTR-mediated apoptosis induced by NGF. Plays a role in zinc-triggered neuronal death. By similarity. May play an important role in the pathogenesis of neurogenetic diseases.

For Research Use Only