

POU6F1 Polyclonal Antibody

Catalog No. E-AB-61631

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human,Mouse,Rat
Immunogen	Recombinant fusion protein of human POU6F1 (NP_002693.3).
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Applications Recommended Dilution

IF	1:50-1:100
-----------	------------

Data



Immunofluorescence analysis of A-549 cells using
POU6F1 Polyclonal Antibody

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

POU6F1 (POU Class 6 Homeobox 1) is a Protein Coding gene. Diseases associated with POU6F1 include Clear Cell Adenocarcinoma Of The Ovary and Clear Cell Adenocarcinoma. Gene Ontology (GO) annotations related to this gene include DNA-binding transcription factor activity. An important paralog of this gene is POU6F2.

For Research Use Only