

TRIAP1 Polyclonal Antibody

Catalog No. E-AB-62447

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human
Immunogen	Recombinant fusion protein of human TRIAP1 (NP_057483.1).
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Applications Recommended Dilution

IF 1:50-1:100

Data



Immunofluorescence analysis of MCF-7 cells using
TRIAP1 Polyclonal Antibody

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

TRIAP1 (TP53 Regulated Inhibitor Of Apoptosis 1) is a Protein Coding gene. Diseases associated with TRIAP1 include Immune-Complex Glomerulonephritis and Autoimmune Glomerulonephritis. Among its related pathways are TP53 Regulates Transcription of Cell Death Genes and Gene Expression. Gene Ontology (GO) annotations related to this gene include p53 binding and phosphatidic acid transporter activity.

For Research Use Only