

CIAO1 Polyclonal Antibody

Catalog No. E-AB-61666

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human,Mouse,Rat
Immunogen	Recombinant fusion protein of human CIAO1 (NP_004795.1).
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.


Applications Recommended Dilution

IHC 1:50-1:200

Data



Immunohistochemistry of paraffin-embedded Rat brain using CIAO1 Polyclonal Antibody at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded Human kidney using CIAO1 Polyclonal Antibody at dilution of 1:100 (40x lens).

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

CIAO1 (Cytosolic Iron-Sulfur Assembly Component 1) is a Protein Coding gene. Diseases associated with CIAO1 include Progressive External Ophthalmoplegia With Mitochondrial Dna Deletions, Autosomal Dominant 6 and Xeroderma Pigmentosum, Variant Type. Among its related pathways are Glucose / Energy Metabolism and Metabolism.

For Research Use Only