

ERGIC2 Polyclonal Antibody

Catalog No. E-AB-61674

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human,Mouse
Immunogen	Recombinant fusion protein of human ERGIC2 (NP_057654.2).
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Applications Recommended Dilution

IHC 1:50-1:200

Data



Immunohistochemistry of paraffin-embedded Mouse kidney using ERGIC2 Polyclonal Antibody at dilution of 1:100 (40x lens).

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

ERGIC2 (ERGIC And Golgi 2) is a Protein Coding gene. Diseases associated with ERGIC2 include Acute Maxillary Sinusitis and Temtamy Syndrome. An important paralog of this gene is ERGIC3.

For Research Use Only