

SPIN2B Polyclonal Antibody

Catalog No. E-AB-62602

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human,Mouse,Rat
Immunogen	Recombinant fusion protein of human SPIN2B (NP_001006682.1).
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Applications Recommended Dilution

IHC 1:50-1:200

Data



Immunohistochemistry of paraffin-embedded Human kidney using SPIN2B Polyclonal Antibody at dilution of 1:100 (40x lens).

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

SPIN2B (Spindlin Family Member 2B) is a Protein Coding gene. Diseases associated with SPIN2B include Rapp-Hodgkin Syndrome and Myasthenic Syndrome, Congenital, 13. An important paralog of this gene is SPIN2A.

For Research Use Only