

OSGEP Polyclonal Antibody

Catalog No. E-AB-61731

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human,Mouse,Rat
Immunogen	Recombinant fusion protein of human OSGEP (NP_060277.1).
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Applications Recommended Dilution

IF 1:50-1:200

Data



Immunofluorescence analysis of A549 cells using
OSGEP Polyclonal Antibody

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

OSGEP (O-Sialoglycoprotein Endopeptidase) is a Protein Coding gene. Diseases associated with OSGEP include Galloway-Mowat Syndrome 3 and Galloway-Mowat Syndrome. Among its related pathways are Gene Expression and tRNA processing. Gene Ontology (GO) annotations related to this gene include transferase activity, transferring acyl groups other than amino-acyl groups and transferase activity, transferring acyl groups. An important paralog of this gene is OSGEPL1.

For Research Use Only