

## ZFAND3 Polyclonal Antibody

Catalog No. E-AB-61734

**Note:** Centrifuge before opening to ensure complete recovery of vial contents.


### Description

<b>Reactivity</b>	Mouse,Rat
<b>Immunogen</b>	Recombinant fusion protein of human ZFAND3
<b>Host</b>	Rabbit
<b>Isotype</b>	IgG
<b>Purification</b>	Affinity purification
<b>Conjugation</b>	Unconjugated
<b>Buffer</b>	PBS with 0.02% sodium azide,50% glycerol,pH7.3.

### Applications Recommended Dilution


**IHC** 1:50-1:200

### Data



Immunohistochemistry of paraffin-embedded rat kidney using ZFAND3 Polyclonal Antibody at dilution of 1:100 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with IHC staining protocol.

**Observed Mw:27kDa**  
**Calculated Mw:25kDa**



Immunohistochemistry of paraffin-embedded mouse kidney using ZFAND3 Polyclonal Antibody at dilution of 1:100 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with IHC staining protocol.

### Preparation & Storage

**Storage** Store at -20°C. Avoid freeze / thaw cycles.

### Background

Predicted to enable DNA binding activity and zinc ion binding activity.

### For Research Use Only