

SMAD3 Polyclonal Antibody

Catalog No. E-AB-62121

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Description


Reactivity	Human, Mouse, Rat
Immunogen	Recombinant protein of human SMAD3
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.

Applications Recommended Dilution

WB 1:500 - 1:2000

IHC 1:50 - 1:200


Data



Western blot analysis of extracts of A-549 cells with SMAD3 Polyclonal Antibody

Observed Mw: 48kDa

Calculated Mw: 25kDa/35kDa/43kDa/48kDa



Immunohistochemistry of paraffin-embedded human kidney with SMAD3 Polyclonal Antibody

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

The protein encoded by this gene belongs to the SMAD, a family of proteins similar to the gene products of the Drosophila gene 'mothers against decapentaplegic' (Mad) and the C. elegans gene Sma. SMAD proteins are signal transducers and transcriptional modulators that mediate multiple signaling pathways. This protein functions as a transcriptional modulator activated by transforming growth factor-beta and is thought to play a role in the regulation of carcinogenesis.

For Research Use Only