

CLPTM1 Polyclonal Antibody

Catalog No. E-AB-62680

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human, Mouse, Rat
Immunogen	Recombinant fusion protein of human CLPTM1 (NP_001285.1).
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Applications Recommended Dilution

WB 1:500-1:2000

Data



Western blot analysis of extracts of various cell lines using CLPTM1 Polyclonal Antibody at dilution of 1:1000.

Observed Mw:100kDa

Calculated Mw:65kDa/75kDa/76kDa

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

CLPTM1 (CLPTM1 Regulator Of GABA Type A Receptor Forward Trafficking) is a Protein Coding gene. Diseases associated with CLPTM1 include Cleft Lip and Cleft Lip/Palate-Ectodermal Dysplasia Syndrome. An important paralog of this gene is CLPTM1L.

For Research Use Only