

GABARAPL1 Polyclonal Antibody

Catalog Number:E-AB-62401



Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human,Mouse,Rat
Immunogen	Recombinant fusion protein of human GABARAPL1 (NP_113600.1).
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Formulation	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Applications Recommended Dilution

WB	1:500-1:2000
-----------	--------------

Data



Western blot analysis of extracts of various cell lines using GABARAPL1 Polyclonal Antibody at dilution of 1:1000.

Observed Mw:16kDa

Calculated Mw:14kDa/17kDa

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

GABARAPL1 (GABA Type A Receptor Associated Protein Like 1) is a Protein Coding gene. Diseases associated with GABARAPL1 include Alzheimer Disease 13 and Vici Syndrome. Among its related pathways are Macroautophagy and Apelin signaling pathway. Gene Ontology (GO) annotations related to this gene include microtubule binding and GABA receptor binding. An important paralog of this gene is GABARAP.

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017