

RASGRP3 Polyclonal Antibody

Catalog No. E-AB-63081

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human,Mouse
Immunogen	Recombinant fusion protein of human RASGRP3 (NP_056191.1).
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Applications Recommended Dilution

IF 1:50-1:200

Data



Immunofluorescence analysis of U2OS cells using RASGRP3 Polyclonal Antibody at dilution of 1:100.
Blue: DAPI for nuclear staining.

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

The protein encoded by this gene is a guanine nucleotide exchange factor that activates the oncogenes HRAS and RAPIA. Defects in this gene have been associated with systemic lupus erythematosus and several cancers.

For Research Use Only