

RAB17 Polyclonal Antibody

Catalog No. E-AB-63165

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human, Mouse, Rat
Immunogen	Recombinant fusion protein of human RAB17 (NP_071894.1).
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.


Applications Recommended Dilution

IHC 1:50-1:200


Data



Immunohistochemistry of paraffin-embedded Rat brain using RAB17 Polyclonal Antibody at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded Human brain cancer using RAB17 Polyclonal Antibody at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded Mouse brain using RAB17 Polyclonal Antibody at dilution of 1:100 (40x lens).

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

RAB17 (RAB17, Member RAS Oncogene Family) is a Protein Coding gene. Diseases associated with RAB17 include

For Research Use Only

Carpenter Syndrome 1. Among its related pathways are Metabolism of proteins and Autophagy Pathway. Gene Ontology (GO) annotations related to this gene include GTP binding and GDP binding. An important paralog of this gene is RAB5A.