

PIP5K1A Polyclonal Antibody

Catalog No. E-AB-62966

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human, Mouse
Immunogen	Recombinant protein of human PIP5K1A
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.

Applications Recommended Dilution

WB 1:500 - 1:2000 IF
1:50 - 1:100


Data



Western blot analysis of extracts of HeLa cells with PIP5K1A Polyclonal Antibody

Observed Mw: 70kDa

Calculated Mw: 56kDa/58kDa/61kDa/62kDa



Immunofluorescence analysis of A549 cells with PIP5K1A Polyclonal Antibody

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

For Research Use Only