

FUBP3 Polyclonal Antibody

Catalog No. E-AB-62292

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human, Mouse
Immunogen	Recombinant fusion protein of human FUBP3 (NP_003925.1).
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Applications Recommended Dilution


WB 1:500-1:2000 IHC
1:50-1:200

Data



Western blot analysis of extracts of various cell lines using FUBP3 Polyclonal Antibody at dilution of 1:1000.

Observed Mw:62kDa
Calculated Mw:28kDa/61kDa



Immunohistochemistry of paraffin-embedded Mouse liver using FUBP3 Polyclonal Antibody at dilution of 1:100 (40x lens).

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

FUBP3 (Far Upstream Element Binding Protein 3) is a Protein Coding gene. Gene Ontology (GO) annotations related to this gene include nucleic acid binding and RNA binding. An important paralog of this gene is FUBP1.

For Research Use Only