

UCHL5 Polyclonal Antibody

Catalog No. E-AB-62446

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human,Mouse,Rat
Immunogen	Recombinant fusion protein of human UCHL5 (NP_001186190.1).
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Applications Recommended Dilution

IF 1:50-1:100

Data



Immunofluorescence analysis of A549 cells using
UCHL5 Polyclonal Antibody

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

UCHL5 (Ubiquitin C-Terminal Hydrolase L5) is a Protein Coding gene. Diseases associated with UCHL5 include Myeloma, Multiple. Among its related pathways are Signaling by GPCR and Metabolism of proteins. Gene Ontology (GO) annotations related to this gene include endopeptidase inhibitor activity. An important paralog of this gene is BAP1.

For Research Use Only