

RUFY2 Polyclonal Antibody

Catalog No. E-AB-63132

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human,Mouse,Rat
Immunogen	Recombinant fusion protein of human RUFY2 (NP_001035882.1).
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Applications Recommended Dilution

IHC 1:50-1:200

Data



Immunohistochemistry of paraffin-embedded Human brain cancer using RUFY2 Polyclonal Antibody at dilution of 1:100 (40x lens).

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

RUFY2 (RUN And FYVE Domain Containing 2) is a Protein Coding gene. Diseases associated with RUFY2 include Dyskeratosis Congenita, Autosomal Recessive 2. An important paralog of this gene is RUFY1.

For Research Use Only