

JMJD4 Polyclonal Antibody

Catalog Number:E-AB-62516

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human,Mouse,Rat
Immunogen	Recombinant fusion protein of human JMJD4 (NP_075383.2).
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Formulation	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Applications Recommended Dilution

WB	1:1000-1:5000
-----------	---------------

Data



Western blot analysis of extracts of various cell lines using JMJD4 Polyclonal Antibody at dilution of 1:5000.

Observed Mw:60kDa
Calculated Mw:50kDa/52kDa

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

JMJD4 (Jumonji Domain Containing 4) is a Protein Coding gene. Diseases associated with JMJD4 include Intracranial Chondrosarcoma and Familial Isolated Hypoparathyroidism. An important paralog of this gene is JMJD6.

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017