

OSGEPL1 Polyclonal Antibody

Catalog No. E-AB-63164

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human, Mouse, Rat
Immunogen	Recombinant fusion protein of human OSGEPL1 (NP_071748.2).
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.


Applications Recommended Dilution

IHC 1:50-1:100 IF
1:50-1:100


Data




Immunohistochemistry of paraffin-embedded Human liver tissue using OSGEPL1 Polyclonal Antibody at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded Rat brain using OSGEPL1 Polyclonal Antibody at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded Mouse stomach using OSGEPL1 Polyclonal Antibody at dilution of 1:100 (40x lens).



Immunofluorescence analysis of A549 cells using OSGEPL1 Polyclonal Antibody

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

For Research Use Only

OSGEPL1 (O-Sialoglycoprotein Endopeptidase Like 1) is a Protein Coding gene. Diseases associated with OSGEPL1 include Galloway-Mowat Syndrome 1 and Mitochondrial Myopathy, Encephalopathy, Lactic Acidosis, And Stroke-Like Episodes. Gene Ontology (GO) annotations related to this gene include metalloendopeptidase activity and transferase activity, transferring acyl groups other than amino-acyl groups. An important paralog of this gene is OSGEP.