

SERPINA12 Polyclonal Antibody

Catalog No. E-AB-63214

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human, Mouse
Immunogen	Recombinant fusion protein of human SERPINA12 (NP_776249.1).
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Applications Recommended Dilution

WB 1:500-1:2000

Data



Western blot analysis of extracts of Mouse brain using SERPINA12 Polyclonal Antibody at dilution of 1:1000.

Observed Mw:47kDa
Calculated Mw:47kDa

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

Adipokine that modulates insulin action by specifically inhibiting its target protease KLK7 in white adipose tissues.

For Research Use Only