

STXBP3 Polyclonal Antibody

Catalog No. E-AB-62935

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human,Mouse,Rat
Immunogen	Recombinant fusion protein of human STXBP3 (NP_009200.2).
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Applications Recommended Dilution

IF 1:50-1:200

Data



Immunofluorescence analysis of U2OS cells using STXBP3 Polyclonal Antibody at dilution of 1:100.
Blue: DAPI for nuclear staining.

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

STXBP3 (Syntaxin Binding Protein 3) is a Protein Coding gene. Diseases associated with STXBP3 include Hemophagocytic Lymphohistiocytosis, Familial, 4. Among its related pathways are Response to elevated platelet cytosolic Ca²⁺ and Translocation of GLUT4 to the plasma membrane. Gene Ontology (GO) annotations related to this gene include syntaxin binding and syntaxin-1 binding. An important paralog of this gene is STXBP1.

For Research Use Only