

IFI44 Polyclonal Antibody

Catalog No. E-AB-63040

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human,Rat
Immunogen	A synthetic peptide of mouse IFI44
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Applications Recommended Dilution


WB	1:500-1:2000
IHC	1:100-1:200
IF	1:50-1:200

Data




Western blot analysis of extracts of Rat thymus using IFI44 Polyclonal Antibody at dilution of 1:1000.


Observed Mw:50kDa
Calculated Mw:47kDa



Immunohistochemistry of paraffin-embedded Human colon using IFI44 Polyclonal Antibody at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded Human stomach using IFI44 Polyclonal Antibody at dilution of 1:100 (40x lens).



Immunofluorescence analysis of U2OS cells using IFI44 Polyclonal Antibody at dilution of 1:100. Blue: DAPI for nuclear staining.

Preparation & Storage

For Research Use Only

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

IFI44 (Interferon Induced Protein 44) is a Protein Coding gene. An important paralog of this gene isIFI44L.

For Research Use Only