

TTLL12 Polyclonal Antibody

Catalog No. E-AB-63069

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human,Mouse,Rat
Immunogen	Recombinant fusion protein of human TTLL12 (NP_055955.1).
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.


Applications Recommended Dilution

IHC 1:50-1:100


Data



Immunohistochemistry of paraffin-embedded Human liver damage using TTLL12 Polyclonal Antibody at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded Rat heart using TTLL12 Polyclonal Antibody at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded Mouse heart using TTLL12 Polyclonal Antibody at dilution of 1:100 (40x lens).

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

TTLL12 (Tubulin Tyrosine Ligase Like 12) is a Protein Coding gene. Diseases associated with TTLL12 include Acute

For Research Use Only

Closed-Angle Glaucoma and Orofaciodigital Syndrome Viii. An important paralog of this gene is TTLL4.