

RSAD2 Polyclonal Antibody

Catalog Number:E-AB-65580

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human,Mouse,Rat
Immunogen	Recombinant fusion protein of human RSAD2 (NP_542388.2).
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Formulation	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Applications Recommended Dilution

IF	1:50-1:200
-----------	------------

Data



Immunofluorescence analysis of U2OS cells using RSAD2 Polyclonal Antibody at dilution of 1:100. Blue: DAPI for nuclear staining.

Preparation & Storage

Storage	Store at -20°C. Avoid freeze / thaw cycles.
----------------	---

Background

RSAD2 (Radical S-Adenosyl Methionine Domain Containing 2) is a Protein Coding gene. Diseases associated with RSAD2 include Yellow Fever and Chikungunya. Among its related pathways are Interferon gamma signaling and Innate Immune System. Gene Ontology (GO) annotations related to this gene include protein self-association and iron-sulfur cluster binding.

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017