

STXBP4 Polyclonal Antibody

Catalog No. E-AB-63225

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Description


Reactivity	Human, Mouse
Immunogen	Recombinant fusion protein of human STXBP4 (NP_848604.3).
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Applications

Recommended Dilution

IHC 1:50-1:100

Data



Immunohistochemistry of paraffin-embedded Human esophagus using STXBP4 Polyclonal Antibody at dilution of 1:100 (40x lens).

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

The STXBP4 antibody from Proteintech is a rabbit polyclonal antibody to a recombinant protein of human STXBP4. This antibody recognizes human, mouse, rat antigen. The STXBP4 antibody has been validated for the following applications: ELISA, WB, IHC analysis.

For Research Use Only