

TAMM41 Polyclonal Antibody

Catalog No. E-AB-62572

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human, Mouse
Immunogen	Recombinant fusion protein of human TAMM41
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Applications Recommended Dilution

WB	1:500-1:2000
-----------	--------------

Data



Western blot analysis of extracts of various cell lines using TAMM41 Polyclonal Antibody at 1:1000 dilution.

Observed Mw: 38 kDa

Calculated Mw: 35 kDa/38 kDa/51 kDa

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

Catalyzes the conversion of phosphatidic acid (PA) to CDP-diacylglycerol (CDP-DAG, an essential intermediate in the synthesis of phosphatidylglycerol, cardiolipin and phosphatidylinositol).

For Research Use Only